



# ROTEX™ 782

## COMPACT ELECTRIC HANDPIECE UNIT

ROTEX™ 782

- Ball bearing high torque electric motor.
- 1,000 to 35,000 R.P.M. variable speed.
- Double insulation heavy duty construction.
- Automatic cooling fan for continuous high speed performance.
- Forward/Reverse switch.
- Standard package includes Variable Speed Foot Control, Handpiece and Handpiece Rest.
- Circuit overload protection.
- Application for grinding, polishing, milling of precious metals, ceramics, acrylics and steel.
- Unit is available for 100/110V 50/60Hz or 220/240V 50Hz operation.

ISO 9001  
ISO 13488



**DENTAMERICA®**

18320 Bedford Circle  
P.O. Box 3200  
Industry, CA 91744 U.S.A.

Phone: (626) 912-1388  
Fax: (626) 913-0510  
www.dentamerica.com

Rev-A-9/16/04