

SCALEX[®] 890

SELF-CONTAINED PIEZO ULTRASONIC SCALER



Self-contained solution tank

SCALEX[®] 890

- Self-contained solution tank allows for convenient portability.
- Any type of irrigant can be used during procedures to meet all clinical needs.
- State of the art piezoelectric self-tuning and load-adjusting technology for smooth and effortless operation.
- Innovative cooling design for greater stability during operation.
- Available selection of a wide variety of high quality tips to correspond with specific clinical applications.
- Unit is available in AC115V/60Hz or 230V/50Hz operation.

ISO 13485

GMP

DENTAMERICA[®]

SCALEX[®] 890

SELF-CONTAINED PIEZO ULTRASONIC SCALER

※PIEZO SCALER TIPS



88002

Reorder No.155-03
Endo Chuck



88003

Reorder No.155-01
Universal A



88004

Reorder No.155-51
Universal B



88005

Reorder No.155-02
Chisel



88006

Reorder No.155-48
Perio



Self-contained solution tank

Easily change the solutions according to different applications including periodontics, endodontics, and prophylaxis.

SPECIFICATIONS:

- FREQUENCY: 27,000 to 32,000 Hz/sec.
- DIMENSIONS: D260xW140xH210 mm
- VOLUME: 1,000 ml
- WEIGHT: 3.0 kg
- INPUT: 115V/60Hz or 230V/50Hz



Handpiece

Autoclavable handpiece with piezo ceramic crystals and titanium tip socket structure.



Foot control

Unique compact foot control designed to be activated at any angle and accessible from 360degrees.

DENTAMERICA[®]

18688 E. San Jose Ave.
City of Industry, CA 91748 U.S.A.
Phone: (626) 912-1388
Fax: (626) 913-0510
www.dentamerica.com

DISTRIBUTOR

DENTAMERICA[®]