

# LITEX<sup>®</sup> 680A

## CURING LIGHT



LITEX<sup>®</sup>  
680A

ISO 13485

GMP

- 75 Watts, 12V, tungsten halogen lamp.
- Internally-mounted super-cool fan for continuous duty.
- Accurately calibrated automatic timer allows preselection of cure rate from 10 to 60 seconds with audible tone.
- Interrupt cycle capability.
- Ergonomically designed “pistol grip” for maximum comfort.
- Autoclavable fiber optic probe, full 360 degree rotation for easy access.
- A thermal switch to safely shut the unit off when overheating.
- Unit is available for 100/110 volts 50/60HTz, or 220/240 volts 50HTz operation.
- Standard package includes 9mm general purpose probe and eye protection shield for safety operation, unit can be wall mounted or free standing.

**DENTAMERICA<sup>®</sup>**

18688 E. San Jose Ave.  
City of Industry, CA 91748 U.S.A.  
Phone: (626) 912-1388  
Fax: (626) 913-0510  
www.dentamerica.com  
www.dentamerica.asia

BC680A  
Rev.C-Mar/11/2022