

LITEX™ 696

CORDLESS LED Curing Light



LITEX™
696

- Thin and lightweight dual intensity light polymerization LED curing light.
- Standard autoclavable 8.5 mm probe with 360° rotation for easy access.
- No more bulbs to replace with average LED life of 50,000 hours.
- Smooth and contoured ergonomic design offers added comfort and easy cleaning.
- Conveniently located thumb buttons set timer and intensity for two stage step-curing.
- Two-stage polymerization with initial lower intensity prior to higher intensity polymerization resulting in less shrinkage stress.
- Automatic timer with audible tone and digital display on handpiece.
- Light intensity of 1,200 mW/cm² cures composite in 10 seconds/4mm.
- Curing cycle can be interrupted manually with one-touch trigger button.
- Large lithium battery capacity provides 750 ten-second curing cycles on a full charge.
- Low battery warning beeper.
- Built-in light intensity indicator.
- Eye protection shield included for safe operation.


ISO 9001
ISO 13485



DENTAMERICA®

LITEX™ 696

CORDLESS LED CURING LIGHT



10-second cure for 4mm curing depth.



Built-in Light Intensity Indicator.



SPECIFICATIONS:

- INPUT: 100-240V AC 0.4A
- OUTPUT: 9V DC 1.3A
- HANDPIECE: 150g
- CHARGING UNIT: 285g
- L180xW36xH40mm
- L134xW92xH157mm

DENTAMERICA®

P.O. Box 3200
Industry, CA 91744 U.S.A.
Phone: (626) 912-1388
Fax: (626) 913-0510
www.dentamerica.com
www.dentamerica.com.tw

DISTRIBUTOR