

SCALEX™ 850

SELF-CONTAINED ULTRASONIC SCALER



Self-contained solution tank

SCALEX™ 850

- Self-contained solution tank allows for convenient portability.
- One touch selection between medicament tank and external water source.
- 25k/30k scaling strokes per second, solid state power unit.
- Innovative water-cooled circuitry designed for greater stability during operation.
- Available selection of a wide variety of high quality inserts to correspond with specific clinical applications.
- Unit is available in 100/110V 50/60Hz or 220/240V 50Hz operation.

ISO 9001
ISO 13485



0434



DENTAMERICA®

SCALEX™ 850

ULTRASONIC SCALER

※ULTRASONIC INSERTS

Reorder No.80005

Universal Insert 25K
(External Water)

**Reorder No.83005**

Universal Insert 30K
(External Water)

**Reorder No.80006**

Universal Insert 25K
(Internal Water)

**Reorder No.83006**

Universal Insert 30K
(Internal Water)

**Reorder No.80007**

Chisel Type Insert 25K
(External Water)

**Reorder No.83007**

Chisel Type Insert 30K
(External Water)

**Reorder No.80008**

Chisel Type Insert 25K
(Internal Water)

**Reorder No.83008**

Chisel Type Insert 30K
(Internal Water)

**Reorder No.80011**

Endodontic Insert 25K
(Canalsonic Tip)

**Reorder No.88060****Foot control**

Unique compact foot control designed to be activated at any angle and accessible from 360° degrees.

SPECIFICATIONS:

- FREQUENCY: 25,000 or 30,000 Hz/sec.
- DIMENSIONS: D260xW140xH210 mm
- VOLUME: 1,000 ml
- WEIGHT: 3.3 kg
- INPUT: 115V/60Hz or 230V/50Hz
- OUTPUT: AC 24V 3A

DENTAMERICA®

P.O. Box 3200
Industry, CA 91744 U.S.A.
Phone: (626) 912-1388
Fax: (626) 913-0510
www.dentamerica.com
www.dentamerica.com.tw

**Reorder No.85003****Self-contained solution tank**

Easily change the solutions according to different applications including periodontics, endodontics, and prophylaxis.

DISTRIBUTOR